




**CARBOFIX  
SPINE**

# CarboClear® Thoracolumbar VBR System

The Perfect Combination of Carbon Fiber & 3D-Printed Endplates



- **Negligible Artifacts in MRI & CT**
  - Enhanced Follow-Up
  - Improved Radiation Planning
- **Optimal Bone In-Growth**
  - Porous Endplates using 3D-Printed Titanium
- **Greater Fatigue** in Comparison to Metal Implants
- **Optimal Radiation Therapy Administration**
  - Negligible Backscattering and Attenuation
- **Extra Large Bone Graft Compartment**
- **Radiolucent**



CarboClear® VBR: MRI



CarboClear® VBR: CT




**CARBOFIX  
SPINE**

# CarboClear® Thoracolumbar VBR System

**Carbon Fiber Implants: Seeing is Believing**

VBR\*

Cat. Number	VBR Type	Cross Section [mm]	Height [mm]	Image
VBR27220XX	Small	22 X 27	30-85 (5mm increments)	
VBR30250XX	Medium	25 X 30	30-85 (5mm increments)	
VBR35275XX**	Large	27.5 X 35	30-85 (5mm increments)	



CarboClear® VBR MRI Images

XX - length, in mm  
\* Supplied Sterile Packed  
\*\* Special Order

**DISTRIBUTED BY:**

CarboFix Spine Inc.  
9983 Traders Lane Ext 1 Calabash,  
NC 28467, USA

**All Implants supplied sterile packed**

Tel: 1-800-408-0120

E-Mail: [info@carbo-fix.com](mailto:info@carbo-fix.com)

**Rx ONLY**

Patents are pending

[www.carbo-fix.com](http://www.carbo-fix.com)